630-012

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UK

INTRODUCTION

To get the best out of your new apparatus, please read this user guide carefully before using it for the first time. Take particular note of the safety precautions. We also recommend that you keep the instructions for future reference.

IMPORTANT SAFETY PRECAUTIONS

- For domestic use only.
- Use for its intended purpose only.
- Do not use if not functioning properly, or if damaged by a heavy impact or similar.
- The armband must not be stretched, twisted or bent.
- Do not use on babies or others who cannot communicate if the armband is too tight, etc.
- Do not immerse in water or any other liquid.
- Do not use in the vicinity of mobile phones or microwave ovens as these emit electro-magnetic radiation that can affect the apparatus.
- Do not expose to direct sunlight, high temperatures, humidity, noise or corrosive substances.
- Never try to repair the apparatus yourself. If defective, return to dealer. Unauthorised repairs or modifications to the apparatus will invalidate the guarantee. Please contact the store where you bought the apparatus for repairs under guarantee.
- This apparatus must not be used to make medical diagnoses. It is suitable only for reference. Consult a doctor if your blood pressure readings seem abnormal.
- Always follow the procedure specified in these user instructions when taking readings.

Batteries

- Remove the batteries if the apparatus is not to be used for a prolonged period.
- Change the batteries immediately once they are dead.
- Never use different types of batteries, or use a new and an old battery at the same time.
- Ensure the batteries are correctly polarised when installing them.

KEY TO MAIN COMPONENTS OF THE APPARATUS

1. Storage box
2. CD
3. USB cable
4. Collar
5. Battery cover
6. Communications jack
7. Display
8. T-button
9. R-button
10. On/off button

Symbols on the display

- a. Hypertension
- b. High normal
- c. Normal
- d. Hypotension
- e. Pulse symbol
- f. Date
- g. Battery symbol
- h. Time/year
- i. Saved reading number
- j. Pulse
- k. Diastolic pressure
- l. Systolic pressure

The date is displayed in the American date display format.

The order is MM-DD-YYYY. For example, 10.12.2006 is to be interpreted as October 12, 2006.

ABOUT BLOOD PRESSURE

What is blood pressure?

When the heart pumps blood around the body, the blood exerts pressure on the walls of the blood vessels. This is known as blood pressure. Blood pressure is different in different parts of the body and depends on the pulse, how elastic the blood vessels are and how thick the blood is. The pressure that can be measured when the heart pumps out into the body is known as systolic blood pressure. The pressure that can be measured when blood runs back from the body towards the heart is known as diastolic blood pressure.

The two types of blood pressure are usually given together, with the systolic blood pressure always given first. If blood pressure is given as "120/70" (expressed as "120 over 70"), this means that the systolic blood pressure is 120, and the diastolic 70.

The reading unit for blood pressure is mmHg (millimetres of mercury).

Why does blood pressure vary?

Blood pressure varies throughout the day and over a person's life. Many different factors affect your blood pressure, something you should be aware of when taking your blood pressure. For instance, be aware of the following:

- Blood pressure is higher during the day than at night.
- Blood pressure increases slightly after you have just eaten. A lot of salt in food can also increase blood pressure.
- Light exercise (e.g. walking and jogging) causes the systolic blood pressure to rise slightly, while hard exercise will cause it to rise dramatically. However, it will fall again when the body relaxes. A lack of exercise raises blood pressure permanently.
- Blood pressure can increase if you are excited, stressed or tense.
- Sudden climatic changes can cause blood pressure to rise.
- The systolic pressure increases with age.
- Women have a greater risk of increased blood pressure after the menopause.

- Insufficient sleep, smoking, obesity, alcohol and certain illnesses cause blood pressure to rise.
- Bathing and toilet visits may cause blood pressure to rise.

When is my blood pressure too high or too low?

High blood pressure is also called hypertension. Low blood pressure is also called hypotension. The diagram below shows when your blood pressure is too high or too low.

Range	Systolic	Diastolic
Hypertension	> 140 mmHg	> 90 mmHg
High normal	130-139 mmHg	85-89 mmHg
Normal	91-129 mmHg	51-84 mmHg
Hypotension	< 90 mmHg	< 50 mmHg

Why should I measure my blood pressure?

Blood pressure naturally increases as you get older, but high blood pressure can also be the first sign of illness, being overweight or of the arteries beginning to stiffen, which increases the risk of blood clots, among other things. It is therefore a good idea to keep an eye on your blood pressure and consult a doctor if it seems unusually high or low.

Important! Do not use to make medical diagnoses. For use only as reference. Consult a doctor if your blood pressure readings seem abnormal. Always follow the procedure specified in these user instructions when taking readings.

PREPARING THE APPARATUS

Inserting the batteries

- Remove the battery cover (5) by pushing it in the direction of the arrow.
- Insert two AAA batteries. Check batteries are correctly polarised (see diagram in battery compartment.) Replace the cover.
- Change the batteries when the battery symbol (g) is shown in the display (7).
- Caution! If the batteries are changed in less than 1 minute, you will not need to reset the clock and any saved blood pressure readings (see below) will remain in the apparatus memory. If it takes more than 1 minute to change the batteries, the display will show all symbols for a moment. The apparatus will emit a short beep and is reset. Any saved blood pressure readings are deleted. The display then shows the reset time.

Setting the clock

- Press and hold the on/off button (10) until the year number flashes on the display. Use the R-button (9) to set year. Press on/off button when setting is correct.
- o If you hold the button down, the display will change faster.
- The month number now flashes on the display. (Please note! The apparatus shows the month before the date). Enter the month number in the same way as the year number, and press the on/off button when it is correctly set.
- The date, hour number and minute number now flash in turn and are set in the same way as described above.
- The display now shows the date (f) and time/year (h).

OPERATING THE APPARATUS

Preparing for a blood pressure reading

- Sit or lie down for 5-10 minutes before the reading.
- Avoid eating, smoking, exercising, bathing or drinking alcohol, coffee or tea before taking the reading.
- Do not go to the toilet immediately before the reading.
- Do not take a blood pressure reading if you are excited, stressed or tense.
- Relax during the reading and do not talk.

Blood pressure reading

- Place the armband (4) on one of your wrists so that the apparatus is on the underside of your arm and tighten it using the Velcro strap.
 - o It does not matter if you take the reading on your left or right wrist.
 - o The armband should be taut but should not pinch.
 - o The side of the armband facing your hand must be 0.5-1 cm from the hand.
- Place the arm to be measured on a pillow on the table with the palm facing upwards and the apparatus at roughly the same height as your heart.
 - o Failure to do so can result in an inaccurate reading.
- Switch on using the on/off button.
- Two beeps will be heard and the number 0 flashes on the display. Shortly afterwards, the armband is pumped up to 185 mmHg, while the pressure is continuously shown on the display.
 - o The apparatus will buzz while the armband is being pumped up.
- When the armband is pumped up to 185 mmHg, the displayed figure falls again and the blood pressure reading starts.
 - o The apparatus emits a number of beeps while the reading is being taken.
 - o The pulse is also measured, and the pulse symbol (e) is shown on the display.
 - o The armband is automatically pumped up if the pressure is too low.
- Once the blood pressure reading has been taken, a long beep is heard and the display shows the systolic pressure (l), the diastolic pressure (k), the pulse (j) and the time.

- There is also a symbol for the range the measured blood pressure falls under: hypotension (d), normal (c), high normal (b) or hypertension (a).
- The results are shown on the display for around 3 minutes and then saved automatically.
 - o Pressing the on/off button will also save. Two beeps will be emitted.
 - o The result will be deleted from the display when saved.
- If you need to take several readings, you should wait 5-10 minutes between each reading.
- o Caution! An individual reading says nothing about your state of health. You can only get an idea of this after a number of readings.

Saving and transferring readings

- The apparatus can save 80 readings. The latest reading is always saved as number 01 in the memory. Older readings are moved down the scale as new ones are saved, and deleted when the memory is full.
- Connect the USB cable supplied (3) to the communications jack (6) to transfer data to a computer. Data can then be shown and processed using the software supplied on CD (2). Refer to the software user guide on the CD.

Displaying saved blood pressure readings

- Press the R-button to show saved readings.
- The first saved reading (01) will show for approx. 4 seconds. The display shows all the data for the reading plus the reading number (i).
- Press the R-button repeatedly to see all subsequent readings. Press and hold the button to scroll through the readings. The reading number is always displayed.
- o You can also press the T-button (8) to quickly scroll through saved readings.
- The apparatus automatically changes back to the time display.

Deleting saved blood pressure readings

- Remove batteries and wait 2 - 3 minutes. Replace batteries. The memory will now be deleted.
 - o Remember to reset the clock.

ERROR MESSAGES

The apparatus registers any problems and shows the following error messages on the display:

Error message	Meaning	Cause	Solution
ER-1	Error during pumping	1. The armband is not correctly attached. 2. The armband is broken. 3. The pump or valves are defective.	1. Attach the armband correctly. 2. Contact the retailer. 3. Contact the retailer.
ER-2	Disturbances or noise during readings reading taken by a doctor.	1. The user moves or speaks during the reading. 2. The user is too nervous.	1. The user must sit still and avoid talking during the reading. 2. Have a blood pressure
ER-3	Error while emptying the armband	1. The valves are defective. 2. The armband is not correctly attached.	1. Contact the retailer. 2. Attach the armband correctly.
ER-4	Abnormal reading	The user may have a heart problem.	Consult a doctor.
HI	The pressure in the armband is too high.	The air tube is blocked, the pump mechanism is defective or the pressure monitor is defective.	Try to take another reading, or contact the retailer.
Battery symbol	Flat batteries	The batteries are almost dead.	Replace the batteries.

STORAGE

When the not in use, store in the accompanying storage box (1). Store in a dry, cool place. Remove the batteries if the apparatus is not to be used for a prolonged period.

CLEANING

When cleaning the apparatus, you should pay attention to the following points:

- Use a soft, dry cloth to clean the apparatus.
- The armband can be cleaned using a damp cloth with a little detergent.
- Do not use strong, corrosive or abrasive cleaning agents for cleaning the apparatus, as these may scratch the surface.

INFORMATION ON DISPOSAL AND RECYCLING OF THIS PRODUCT

Please note that this Adexi product is marked with this symbol: 

This means that this product must not be disposed of together with ordinary household waste, as electrical and electronic waste must be disposed of separately.

In accordance with the WEEE directive, every member state must ensure correct collection, recovery, handling and recycling of electrical and electronic waste. Private households in the EU can take used equipment to special recycling stations free of charge. In some member states you can, in certain cases, return the used equipment to the retailer from whom you purchased it, if you are purchasing new equipment. Contact your retailer, distributor or the municipal authorities for further information on what you should do with electrical and electronic waste.

GUARANTEE TERMS

The guarantee does not apply:

- if the above instructions are not followed
- if the apparatus has been interfered with
- if the apparatus has been mishandled, subjected to rough treatment, or has suffered any other form of damage

Due to the constant development of our products in terms of function and design, we reserve the right to make changes to the product without prior warning.

IMPORTER

Adexi Group

www.adexi.eu

We cannot be held responsible for any printing errors.

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